Aim: Subjects: To show how the voltage of a charged capacitor changes when changing the dielectric. 5C20 (Dielectric)

5C30 (Energy stored in a capacitor)

**Diagram:** 



#### **Equipment:**

Parallel plate capacitor.

- Piece of isolating foam between the plate capacitor and table.
- Electrostatic voltmeter, 0-25kV.
- Glass plate, thickness 20mm.
- Perspex container,  $50 \times 30 \times 3 cm^3$ , filled with <u>demineralised</u> water.
- Wooden bar as a guiding construction (very well visible in Figure 2).
- Power supply (Leybold Didactic 52237), 0-6kV (see Safety). For the demonstration it is better to use a 25kV power supply, set at 15kV
- Use connection leads with Teflon isolation.
- Electric heater (to prevent that moisture spoils the demonstration).

#### Safety:

• The 6kV power supply has passive current limitation, ensuring that no dangerous contact voltage can occur. Nevertheless be careful, because when you accidentally touch the high, you probably make a spastic movement, and you may hurt yourself (or somebody else).

<sup>\*</sup> In accordance with IEC 61010-1 (Safety requirements for electrical equipment), a part is not deemed to be live (i.e. carrying a dangerous contact voltage) when, at voltages greater than extra-low voltage (> 60 V DC), the current through an induction-free resistance of 2 kOhm is not greater than 2 mA for DC, additionally, the charge for voltages up to 15 kV is less than 45  $\mu$ C, and the stored energy does not exceed 350 mJ for voltages over 15 kV.



**Presentation:** The setup of the demonstration is explained to the students. The plates are set at such a separation that *d* just a little bit larger than the thickness of the glass–plate. The power supply is set at 15kV.



The capacitor is charged by shortly touching the capacitor with the free lead of the 15kV power supply (see Figure 1). After this charging of the capacitor, the voltmeter reads 15kV.

The students are asked what will happen to the voltage of a charged capacitor when the glass plate is shifted between the plates (see Figure 2). After their answers shift the glass between the capacitor plates. They will see that the voltage lowers. Shift carefully all the time, with the glass plate sliding along the grounded plate of the capacitor, so the glass plate does <u>not touch the high voltage positive plate!</u>



Figure 2 Removing the glass plate will increase the voltage again to its original value of 15kV.



The same experiment is performed with the container filled with water (see Figure 3).



**Explanation:** <u>Explanation 1:</u> The demonstration shows that the voltage lowers when a dielectric is shifted between the capacitor plates. Voltage is the energy per unit charge, so we can say that the potential energy of the capacitor lowers. Where is this energy gone? It was needed to polarize the dielectric! The "lost" energy is now in the polarized dipoles of the dielectric.

In the same way the energy stored in the capacitor becomes higher when the dielectric is removed. Also the demonstrator has to pull now because there is an attracting force between the induced charge on the dielectric and the plates. So he adds energy to the capacitor.

Explanation 2: When the dielectric is between the plates, C increases ( $C = \varepsilon \frac{A}{d}$ )

because  $\varepsilon$  becomes higher, and since Q is constant this will cause V to decrease (Q=CV).

(This explanation is just mathematics and no physics.)

Especially when water is the dielectric the voltage should lower dramatically, since  $\mathcal{E}_r$  of water is around 80!

**Remarks:** 

• In the explanation of the demonstration it is supposed that the charge on the capacitor is constant. This is true only if the capacitance of the electrostatic voltmeter is small compared to that of the parallel plate capacitor.



Sources:

<u>Mansfield, M and O'Sullivan, C., Understanding physics</u>, pag. 455-458.
<u>Giancoli, D.G., Physics for scientists and engineers with modern physics</u>, pag. 621-625.

