


Resolution

Aim: To show how diffraction limits the resolution of an optical system.

Subjects: 6C10 (Diffraction From Two Sources)

Diagram:



Equipment:

- Rotatable disc with 8 holes: 3.0, 2.5, 2.0, 1.5, 1.0, 0.5, 0.4 and 0.3 mm (see Figure1).
- Aluminium foil with 2 pair of holes fitted on a stand (see Diagram and Figure2).
- Lamp, 220V/200W.
- Variable transformer on the 220V line voltage.
- Camera with zoom lens.

Safety:

- No remarks.

Resolution

Presentation: Preparation


Built the demonstration as shown in the Diagram:

- Focus the camera on the Aluminium foil.
- The rotatable disc is placed as close as possible to the camera.
- The lamp should not be too close to the Aluminium foil, because we need parallel light beams from the holes in the Aluminium foil. To avoid scattered light a cardboard tube is placed between lamp and the Aluminium foil.
- Adjust the vertical and horizontal position of the lamp and also its intensity to get a satisfying illumination of the small holes in the foil.



Figure 1

The lamp is switched on. The rotatable disc has its largest hole in position. The camera is focussed at the pairs of holes in the aluminium foil. The holes of both pairs in the aluminium foil are observed as separate images.



diameter of holes:
~ 0.5-1.0mm

Figure 2

Presentation

Choose the smallest hole in the rotatable disc. Two rather hazy patches of light are observed by the camera. One of the two patches gives the idea that it could be a double spot. Then a larger hole is selected on the rotatable disc and we see that our idea of one of the light patches being two separate spots is strengthened.

When we continue to select larger holes on the rotatable disc the light spot resolves as really consisting of two light spots. Even the other light spot finally resolves into two! Figure 3 shows the sequence of the observed light spots.

Resolution






					
Diaphragm (mm)	.3	.4	.5	1.0	1.5
$\frac{1.22\lambda}{D}$.002	.0015	.0012	.00061	.00041

Figure 3

(In demonstrating we also go again backwards to smaller holes in the diaphragm.)

Explanation: If two point objects are very close, the diffraction patterns of their images will overlap. As the objects are moved closer, a separation is reached where you can't tell if there are two overlapping images or a single image. The separation at which this happens is stated by Lord Rayleigh: *two images are just resolvable when the centre of the diffraction disk of one image is directly over the first minimum in the diffraction disk of the other*. A circular hole shows a diffraction pattern with a central maximum of half width: $\theta = \frac{1.22\lambda}{D}$, where D is the diameter of the circular opening. Calculating with $\lambda = 500\text{nm}$ we get for the smallest hole on the rotatable disc $D = .3\text{mm}$, $\theta = 2 \times 10^{-3}$.




Figure 4

Resolution

In our demonstration: $\theta = \frac{\Delta l}{f}$ (see Figure 4).

The distance $f = 2$ meter, and calculating θ for both pair of holes, we get $\theta = .75 \times 10^{-3}$ for the pair of holes with a separation of 1.5 mm and $\theta = 2 \times 10^{-3}$ for the pair of holes with a separation of 4 mm.

These calculations compared with the Rayleigh criterion (that is expressed as $\theta = \frac{1.22\lambda}{D}$ and is calculated and listed in the bottom row of Figure 3), shows that the two holes with a separation of 4 mm will be resolved when the diaphragm is larger than .3 mm and that the holes with a separation of 1.5 mm will be resolved when the diaphragm is larger than .5 mm. The observed light spots in Figure 3 show that this is more or less right!

Remarks:

Since $\theta = \frac{1.22\lambda}{D}$, it is useful to do this demonstration in different colours. (We didn't try this yet.)

Sources:

- [Giancoli, D.G., Physics for scientists and engineers with modern physics](#), pag. 896-899
- [Hecht, Eugene, Optics](#), pag. 416 and 461-465
- [PSSC, College Physics](#), pag. Laboratory Guide, Experiment 22