


Water dropper

Aim: To show the principle of operation of a Van der Graaff generator.

Subjects: 5A40 (Induced Charge)
5A50 (Electrostatic Machines)

Diagram:



- Equipment:
- Electrostatic Voltmeter
 - Metal beaker on top of electrostatic voltmeter.
 - Metal collar, fitted on a bar of PVC.
 - Metal container with pinchcock, allowing water to drip. (The holding bar is partly made of perspex for isolation.)
 - Power supply (about 1kV).
 - Camera.
 - Electrostatic Voltmeter (see Remarks and Figure2).

Water dropper

Presentation: Set up the equipment as shown in the Diagram (the metal collar is in the beginning turned away from the dropper; the +1kV lead is not yet connected). The power supply is adjusted to a voltage of about 1kV. Touching momentarily the electroscope with the +1kV lead shows that this voltage gives a small deflection (in our case 3 divisions). Discharge the electroscope.

Now the upper container is filled with water and placed about 10 cm above the beaker on the electroscope. The +1kV lead is connected to it. The pinchcock is opened and drops of water fall into the beaker on the electroscope. When the upper container is connected to the power supply, the electroscope will build up the small deflection (of 3 divisions) we had before.

Now the metal collar, connected to ground, is placed around the falling drops of water. The effect is that the deflection of the electroscope increases and very soon runs off scale!

When now the collar is moved away from the falling drops of water, the electroscope will slowly return to the lower deflection it had before.

Explanation: In the first part of the demonstration an electric field (1kV over 10cm, so $E=10\text{kV/m}$) exists between the upper watercontainer and the electroscope.




Figure 1

Due to this field, the drops of water will be charged when they break loose from the metal dropper (see Figure1). When the drops fall into the beaker on the electroscope this charge accumulates (on the outside of the beaker) and the deflection of the electroscope increases. This continues until the potential of the beaker is the same as the potential of the upper watercontainer, because now there is no longer an electric field to charge the drops of water.

When the grounded metal collar is placed around the stream of falling drops, there will be an electric field again between the collar and the upper watercontainer and again the drops of water will be charged, and again the deflection of the electroscope increases. As long as drops of water fall, the charge on the electroscope increases!

When now the collar is removed, the electric field between the beaker and container changes direction. So, a drop of water will be charged opposite to the charge it got in the first part of the demonstration, lowering the deflection of the electroscope until the electric field between beaker and container is zero again.

Water dropper

Remarks:

- By means of a camera the (reading of the) electrostatic Voltmeter is projected on a monitor or large screen.
- The electrostatic Voltmeter can be replaced by an electrostatic Voltmeter (see Figure2). Then the students can see directly that the voltage of the metal beaker can reach a voltage much higher than the voltage of the container from where the drops are falling.




Figure 2